

Who's Guarding You & Your Animals?

Peroxigard® Ready-to-Use and Wipes One-Step Disinfectant Cleaner for use in Life Sciences

One technology for all your facility's cleaning and disinfection needs with no dilution required.

Efficacy

Effective against over 90+ key pathogens, including those on bioexclusion lists. Achieves disinfection in only 1 minute* against bacteria, viruses, TB, and fungi.

Safety

Eliminates the need to interact with legacy hazardous concentrates. Designed with animal safety in mind. Fragrance free, non-irritating to eyes and skin, non-corrosive and superior compatibility. Leaves no active residue on surfaces or equipment.

Simplicity

Removes the risks associated with diluting concentrates and eliminates the complexity of using multiple products at various dilution ratios. No mixing required and a long 3 year shelf life, reducing your staff's time.



Disinfection Compliance for Better Science

Peroxigard Ready-to-Use (RTU) and Wipes utilize Accelerated Hydrogen Peroxide®, a patented and synergistic blend of commonly used ingredients, that when combined with low levels of hydrogen peroxide create an ideal one-step cleaner and disinfectant. Trusted EPA Registered disinfectants with proven efficacy claims and excellent safety profile.

Peroxigard RTU and Wipes provide an easy, convenient and effective way to disinfect your surfaces and prevent cross contamination. Ideal for work stations, lab benches, electronic equipment (screens, phones), and other small surfaces. Suitable for non-human primate facilities with tuberculocidal activity in only 1 minute. Peroxigard RTU and Wipes are proven to remain effective in the presence of a 5% soil load, providing disinfection confidence.

Effective against key pathogens in 1 minute* or less including:

Bacteria

Klebsiella pneumoniae; *Bordetella bronchiseptica*;
Methicillin-resistant Staphylococcus aureus (MRSA);
Salmonella enterica; *Corynebacterium bovis* (C. bovis)

Viruses

Murine Parvovirus; *Murine Norovirus*; *HIV Type 1*; *Rotavirus*;
Murine Hepatitis Virus; *Adenovirus*; *SARS-CoV-2* (Covid-19 Virus)

Tuberculocidal

Mycobacterium bovis

Fungi

Candida albicans; *Microsporium canis*; *Trichophyton interdigitale*



FOR USE IN:

- Vivariums
- Breeding facilities
- Lab animal research facilities
- Laboratories
- Non-human primate facilities

FOR USE ON:

- Biosafety cabinets
- Lab benches and exam tables
- Diagnostic and laboratory equipment
- Cages, crates, and carts

SHELF LIFE

- 3 years (Opened/Unopened)

AVAILABLE FORMATS

- Ready-To-Use Liquid (EPA Reg No. 74559-9)
32 oz, 1 gal, 5 gal, 30 gal, 55 gal
- Presaturated Wipes (EPA Reg No. 74559-10)
160 ct.

To learn more about Accelerated Hydrogen Peroxide or check out our peer-reviewed studies, please visit: Peroxigard.com



Peroxigard® is a member of the Virox family of brands.

DISCLAIMER: You must read and follow the instructions on the label for our product and use the product for the particular applications specified on the product label. To the full extent permitted under applicable law, Virox Technologies Inc. shall not be liable for any damages whatsoever, arising out of or otherwise in connection with any failure to follow those instructions or any use of its products which is not specified on the product's label. The instructions on the product label are not a warranty, express or implied, and Virox Technologies, Inc. expressly disclaims any implied warranties, including any warranty of fitness for a particular purpose. Under applicable law, some of the above waivers or disclaimers of warranty, or exclusions of liability for certain types of damages may not apply to you. In such cases, our liability will be limited to the fullest extent permitted by applicable law. Product shots have been digitally altered and are not accurate as shown. Refer to product label or reference sheet for a complete list of disinfection claims and additional features. *See reference sheet for a complete list of efficacy claims, additional features and directions for use. Use a 2-minute contact time for *C. bovis*. AHP® and Design and Accelerated Hydrogen Peroxide® and Design are trademarks of Diversey, Inc.